

## RNAtag-Seq Experimental Protocol

# Simultaneous generation of many RNA-seq libraries in a single reaction

March 2015

Please cite: Shishkin et al. *Nature Methods* 2015 2015 Apr;12(4):323-5. doi:  
10.1038/nmeth.3313

Updated versions of this protocol can be found by visiting the Guttman lab website at  
[www.lncRNA.caltech.edu](http://www.lncRNA.caltech.edu)

## RNAtag-seq Protocol for batch of 32 RNA samples

### Reagents:

- Turbo DNase; Ambion/Applied Biosystems, Cat.# AM2239
- FastAP Thermosensitive Alkaline Phosphatase; Thermo Scientific, Cat.# EF0651
- RLT buffer (RNeasy Lysis Buffer); Qiagen, Cat.# 79216 (220mL)
- T4 RNA Ligase 1 (Custom Order from NEB, contact Kellie Dodd); 30,000U/mL (0.5mL) - 3x concentrated, normally 10,000U/mL, Cat.# MO 204L and comes with the following:
  - ATP (100mM); alternative vendor: Roche, Cat.# 11140965001 (aliquot and store at -80°C, always use fresh aliquot)
  - PEG 8000, 50% in water
- DMSO (100%); Sigma, Cat.# D8418-50ML for Molecular Biology
- RNase Inhibitor, Murine (40U/μL); NEB, Cat.# M0314S (3,000 units)
- SUPERase-IN 20U/μL (Ambion; AM2694; 2500U)
- AffinityScript Multiple Temperature cDNA Synthesis Kit, 50 reaction; Agilent, Cat.# 200436 – includes the dNTPs, 10x RT buffer, RNase Block Ribonuclease Inhibitor (40U/μL)
- 10X AffinityScript RT buffer; Agilent, Cat.# 600100-52
- RNA Clean & Concentrator™-5 columns; Zymo Research, Cat.# R1015 (50 preps)
- RNAClean XP beads; Agencourt/Beckman, Cat.# A63987 (40mL)
- Ribo-Zero:
  - Bacteria:  
Ribo-Zero™ Magnetic Gold Kit (Bacteria); Epicentre/Illumina, Cat.# MRZMB126
  - Human:  
Ribo-Zero™ Magnetic Gold Kit (Human/Mouse/Rat); Epicentre/Illumina, Cat.# MRZG126
- Sodium hydroxide solution (5M); Sigma-Aldrich, Cat.# 656046 Fluka
- Glacial acetic acid (17.4M); Sigma-Aldrich, Cat.# ARK2183
- Nuclease free water; Ambion, Cat.# 4387936
- 1x low TE (10mM Tris pH 8.0, 0.1mM EDTA); Affymetrix/USB Cat.#75793
- 2100 Electrophoresis Bioanalyzer Instrument; Agilent, Cat.# G2939AA
- Agilent RNA 6000 Pico Kit; Agilent, Cat.# 5067-1513
- QIAvac 24 Plus Vacuum manifold (1-24 spin columns; Qiagen cat# 19413)

### Plates/Tubes/Caps:

- Flat PCR 12-cap strips, Optically clear; USA Scientific, Cat.# 1400-3120
- TempPlate No-skirt PCR Plates, 96 wells, 10 plates; USA Scientific, Cat.# 1402-9608
- MaxyClear; Maxymum Recovery; 0.2mL; Flat cap (Axygen, cat#PCR-02-L-C)

### Procedure:

#### **1. Quality and quantity of RNA with Agilent Bioanalyzer**

- Check RNA quality by running on the Agilent Bioanalyzer
- Calculate up to 400ng and place in a tube.
- Increase the volume to 15μL with Nuclease free water
- Add 1uL of SUPERase-IN (20U/μL)

- Final total volume = 16 $\mu$ L (25ng/ $\mu$ L)
- Continue to next step or freeze at -80 until ready to process samples

## 2. Fragmentation and DNase –FastAP combination step

### 2.1.Fragmentation using 2x FastAP buffer

- Set up the following reaction:

Reagents	1 reaction	Final conc
FastAP buffer (10x)	4 $\mu$ L	2x
RNA (from step 1.)	16 $\mu$ L	up to 500ng total RNA
Total vol ( $\mu$ L)	20 $\mu$ L	

- Add 4 $\mu$ L of FastAP buffer to the RNA and mix well.
- Incubate on pre-heated thermal cycler for **3 min at 94°C**
  - a. **If RNA is partially degraded (RIN<7), fragment 3min at 92°C – prevents over-fragmenting samples**
- Place on ice

### 2.2.DNase and FastAP treatment

- Set up the following reaction mix

Reagents	1 reaction	Final
fragmented RNA (from step 2.1)	20 $\mu$ L	1x
RNase Inhibitor, Murine (40U/ $\mu$ L)	1 $\mu$ L	10U
TURBO DNase (2U/ $\mu$ L)	4 $\mu$ L	8 U
FastAP (1U/ $\mu$ L)	10 $\mu$ L	10U
Nuclease free water	5 $\mu$ L	
Total vol ( $\mu$ L)	40 $\mu$ L	

- Make enzyme/water master mix and mix well
- Aliquot 20  $\mu$ L into each tube/well with the fragmented RNA
- Mix well
- Incubate on preheated thermal cycler for 30 min at 37°C

### 2.3.SPRI (2x)

- Increase reaction volume to 80  $\mu$ L with nuclease free water (**add 40  $\mu$ L water/reaction**)
- Add 2.0x reaction volume of Agencourt RNAClean XP beads (**160  $\mu$ L**) and follow manufacturer's protocol to capture RNA on beads: Briefly
  - Incubate at room temperature for 15min to bind RNA
  - Place on magnet for 5 min or until solution is clear
  - Aspirate out and discard clear solution
  - Add 200  $\mu$ L fresh 70% EtOH without removing from magnet and incubate for 30sec
  - Aspirate out and discard the clear solution
  - Add 200  $\mu$ L fresh 70% EtOH without removing from magnet and incubate for 30sec
  - Aspirate out and discard the clear solution and let air dry for 10min
- Elute off beads with **12  $\mu$ L nuclease free water**
- Take 5  $\mu$ L of each sample and proceed to 1<sup>st</sup> ligation
- **QC: Save 1.2  $\mu$ L from remaining RNA before addition of SUPERase-IN**
  - Run 11 random samplings on Agilent RNA pico chip to check the fragmentation profile of each batch of 32

- Add 1 $\mu$ L SUPERase-IN (20U/ $\mu$ L) to the remaining material and store at -80°C.
  - this can be used as back up if it is necessary to repeat process

### 3. FIRST LIGATION (RNA/DNA) 3' Linker Ligation

**NOTE:** Set up RNA/barcoded adapter in single tube or use the TempPlate No-skirt PCR Plates for batched samples (these rigid plates are easy to handle and can mix well by flicking by hand) – with a razor, cut a column (for 8 samples) or row (for 12 samples) and use as strip and cover with the Flat PCR 12-cap strips (these strips fit tightly on these plates and will not leak)

- Add 1  $\mu$ L of barcoded adapter (100 pmole\*) to 5  $\mu$ L of dephosphorylated RNA (\*100 pmole = 1  $\mu$ L of 100 $\mu$ M; **barcoded adapters need: 5'-P and 3' ddC (or 3'-C3 spacer)**)
- **Heat at 70°C for 2 min** →place **in cold block** on ice
- Set up Ligation mix below

**NOTE:**

- Make up mix at **room temp** so the reagents don't start precipitating when combined (if DMSO is added directly into cold buffer it will precipitate)
- Pipette very slowly for accurate aspiration of PEG (very viscous)
- When setting up mix for multiple reactions include 25% extra to account for pipetting error due to the viscosity of the PEG in the mix
- All reagents except enzymes (-20°C ) should be stored at -80°C in single use aliquots:

Ligation Mix (set up at room temp)	1 reaction	40 reaction mix (for set of 32)
10 $\times$ T4 RNA Ligase Buffer	2 $\mu$ L	80 $\mu$ L
DMSO (100%)	1.8 $\mu$ L	72 $\mu$ L
ATP (100 mM)	0.2 $\mu$ L	8 $\mu$ L
PEG 8000 (50%)	8 $\mu$ L	320 $\mu$ L
20-50ng of Repaired RNA (tested up to 250ng)	5	
Barcoded Adapter, 100 pmoles (1 $\mu$ L of 100 $\mu$ M stock)	1 (5 $\mu$ M final)	
RNase inhibitor, Murine (40U/ $\mu$ L)	0.3 $\mu$ L	12 $\mu$ L
T4 RNA Ligase 1 (30,000U/mL)	1.8 $\mu$ L	72 $\mu$ L
Total	20.1 $\mu$ L	564 $\mu$ L

- **Mix really well by flicking tube since the solution is very viscous,**
- **Quick spin of mastermix tube**
- Add 14.1  $\mu$ L of ligation mix to each well containing 6  $\mu$ L denatured RNA + adapter
- **Mix well MANY times; mix by flicking since the solution is very viscous**  
(If setting this up on a robot, the reaction should be mixed up/down for 5-10min with low retention tips to mix well!)
- Incubate at 22°C (room temp) for 1 hour 30 minutes

### 4. Pooling step: RLT buffer + Zymo column

- Maximum binding capacity of columns is 5  $\mu$ g; do not exceed when pooling samples
- Add 60  $\mu$ L of **RLT buffer** to each sample to inhibit ligase and mix well (80  $\mu$ L total)
- Pooling and clean up using Zymo Clean & Concentrator™-5 column (200nt cut off)
  - Follow manufacturer's 200nt cut off protocol
  - Add 2x reaction vol (160  $\mu$ L=2x 80  $\mu$ L) of 1:1 binding buffer: EtOH (100%) / reaction

**NOTE:** When pooling multiple samples (>700  $\mu\text{L}$ ) onto Zymo column use a vacuum manifold then proceed to centrifugation steps according to the manual

- Elute 2 times with 16  $\mu\text{L}$  nuclease free water for a total volume of 32  $\mu\text{L}$
- **NOTE:** 2 elutions help improve recovery/yield of RNA
- **Optional: Save 2  $\mu\text{L}$  for QC-** Run on Agilent RNA pico chip to check the fragmentation profile

## 5. Quick Protocol for Ribo-Zero™ See Appendix (bacterial kit)

### 6. First Strand cDNA

- Take **12  $\mu\text{L}$**  rRNA depleted RNA (use all the material from above )
- Add **2  $\mu\text{L}$**  (50 pmoles) of AR2 primer (25  $\mu\text{M}$ )
  - 5' TAC ACG ACG CTC TTC CGA T 3' — AR2 – 53% GC, 19bp,
- Mix well
- **Heat the mixture to 70°C for 2 min and immediately place on ice** (on cold block on ice)
- Make RT mix below: For multiple samples, a mastermix can be prepared ahead of time and added to the RNA/AR2 tube on ice.

Add (in order on ice):

RT Mix	1 reaction
10× AffinityScript RT Buffer	2 $\mu\text{L}$
DTT (0.1M)	2 $\mu\text{L}$
25mM dNTP Mix (25mM each)	0.8 $\mu\text{L}$
RNase inhibitor, murine (40U/ $\mu\text{L}$ )	0.4 $\mu\text{L}$
AffinityScript RT Enzyme	0.8 $\mu\text{L}$
rRNA depleted RNA + AR2 primer	14 $\mu\text{L}$
Total	20 $\mu\text{L}$

- Add 6  $\mu\text{L}$  of RT mix to the 14  $\mu\text{L}$  rRNA depleted RNA + AR2 RT primer on ice
- Mix well and spin for 5 sec
- Place in **HOT** (55°C) incubator or thermocycler. **Incubate at 55°C for 55 minutes**

### 7. RNA degradation after Reverse Transcription

**NOTE:** make fresh working stock solutions of NaOH and Acetic Acid

- Add 10% reaction vol. of 1N NaOH (**2  $\mu\text{L}$** ) to each reaction
- Incubate at 70°C for 12 minutes
- Neutralize with 4  $\mu\text{L}$  of 0.5M Acetic Acid; **mix well**
- Total volume = 26  $\mu\text{L}$

### 8. Reverse Transcription primer cleanup (2x SPRI)

- Add 14  $\mu\text{L}$  of sterile water to each reaction for a final volume of 40  $\mu\text{L}$
- **Transfer to new tubes** (NaOH will start degrading tubes)
- Add 2x reaction volume RNAClean XP beads (80 $\mu\text{L}$ ) to the sample in new tubes, and mix up/down 15x

- Incubate at room temperature for 15min
- Place on magnet for 5 min or until solution is clear
- Aspirate out and discard clear solution
- Add 200  $\mu\text{L}$  fresh 70% EtOH without removing from magnet and incubate for 30sec
- Aspirate out and discard the clear solution
- Add 200  $\mu\text{L}$  fresh 70% EtOH without removing from magnet and incubate for 30sec
- Aspirate out and discard the clear solution and let air dry for 10min
- Elute off beads with 5  $\mu\text{L}$  RNase/DNase free water - **DO NOT TRANSFER TO NEW TUBES AND KEEP BEADS**

### 9. Second LIGATION (ssDNA/ssDNA) 3' Linker Ligation with beads

-3Tr3 adapter: 5' AGATCGGAAGAGCACACGTCTG 3' = 55% GC, 22bp; Needs: 5'-P and 3' ddC (or 3'-C3 spacer)

- Add 2  $\mu\text{L}$  of 3Tr3 adapter to cDNA
- Heat at 75°C for 3 min; Place on cold block on ice
- Make ligation reaction mix: For multiple samples, a mastermix can be prepared ahead of time and added to the cDNA/3Tr3 tube on ice.

2 <sup>nd</sup> Ligation Mix	1 reaction
10 $\times$ T4 Ligase Buffer	2 $\mu\text{L}$
DMSO (100%)	0.8 $\mu\text{L}$
ATP (100 mM)	0.2 $\mu\text{L}$
PEG 8000 (50%)	8.5 $\mu\text{L}$
cDNA	5 $\mu\text{L}$
3Tr3 adapter, 80 pmoles (2 $\mu\text{L}$ of 40 $\mu\text{M}$ stock)	2 $\mu\text{L}$ (4 $\mu\text{M}$ final)
T4 RNA Ligase 1 (30,000U/mL)	1.5 $\mu\text{L}$
Total	20 $\mu\text{L}$

- Swirl the cDNA/beads/water with pipet tip PRIOR to dispensing 13  $\mu\text{L}$  ligation mix
- Mix well by pipetting up and down  $\sim 10\text{x}$  or cap tubes and flick several times; solution is viscous
- Quick spin (low speed centrifuge, to get everything to bottom of tube)
- Incubate overnight at 22°C

### 10. SPRI clean up (2x): Add fresh SPRI beads

- Increase to 2x reaction volume with water (to dilute out viscous ligation reaction)
- Add 2x reaction vol. (80  $\mu\text{L}$ ) of fresh SPRI beads
- Mix up/down 15x
- Let sit at room temperature for 15 minutes
- Place on magnet 5 min or until solution is clear and aspirate and discard clear solution
- Add 200  $\mu\text{L}$  fresh 70% EtOH without removing from magnet and incubate for 30 sec
- Aspirate out and discard the clear solution
- Add 200  $\mu\text{L}$  fresh 70% EtOH without removing from magnet and incubate for 30 sec
- Aspirate out and discard the clear solution
- Let air dry for 10min

- Elute off beads by pipetting 25  $\mu$ L RNase/DNase free water and transfer to new tube.

### 11. 2<sup>nd</sup> SPRI clean up (1.5x) to remove the remaining adapters

- Add 1.5x reaction volume RNAClean XP beads  
25  $\mu$ L x 1.5 = 37.5  $\mu$ L SPRI /reaction
- Mix up/down 15x
- Let sit at room temperature for 15 minutes
- Place on magnet 5min (or until solution is clear). Aspirate out and discard clear solution
- Add 200  $\mu$ L fresh 70% EtOH without removing from magnet and incubate for 30 sec
- Aspirate out and discard the clear solution
- Add 200  $\mu$ L fresh 70% EtOH without removing from magnet and incubate for 30 sec
- Aspirate out and discard the clear solution
- Let air dry for 2-3 min
- Elute with 25  $\mu$ L RNase/DNase free water

### 12. PCR Enrichment-TEST to determine final cycle number

- Set up a test PCR using 5  $\mu$ L of ss cDNA sample and 9-12 cycles of PCR (based on experience with pool of 16 reactions, each starting with ~400ng total RNA)
- Include a negative control for each primer set

- Make a mix consisting of:

PCR Mix	1 reaction
water	15.4 $\mu$ L
10x AccuPrime PCR buffer 1	2.5 $\mu$ L
Primer 1 (2P_univP5, 12.5 $\mu$ M)	1 $\mu$ L
Primer 2 (2P_barcode, 12.5 $\mu$ M)	1 $\mu$ L
cDNA	5 $\mu$ L
AccuPrime HiFi Taq (5U/ $\mu$ L)	0.1 $\mu$ L
<b>Total</b>	25 $\mu$ L

- Mix well
- Aliquot 18  $\mu$ L / well
- Add 2  $\mu$ L total of primer to each well
- Add 5  $\mu$ L of ss cDNA or water (for negative control)
- Mix well and aliquot 8ul into 1 well of the 384 well plate(s) x 3 plates
- Place each plate in a thermal cycler with the following conditions

Cycling conditions:

1 cycle: 98°C, 3min

cycles: 9, 12, 15 cycles (for test PCR)

98°C, 30sec

55°C, 30sec

65°C, 60sec

1cycle: 65°C, 10min

4°C, hold

### 13. SPRI clean up (1.5x):

- Increase to 25uL volume with sterile water
- Add 1.5x AMPure beads (37.5uL), and mix up/down 15x
- Incubate at RT for 15min
- Place on magnet for 5min or until solution is clear

- Aspirate out and discard clear solution
- Add 200  $\mu\text{L}$  fresh 70% EtOH without removing from magnet and incubate for 30 sec.
- Aspirate out and discard the clear solution
- Add 200  $\mu\text{L}$  fresh 70% EtOH without removing from magnet and incubate for 30 sec
- Aspirate out and discard the clear solution and let air dry for 10min
- Elute off beads with 10  $\mu\text{L}$  1x low TE (10 mM Tris, 0.1M EDTA)

#### 14. QC test enrichments on Agilent DNA HS chip

-Based on test results change the cycle number, if necessary, and set up more reactions to provide enough material to send for sequencing

#### 15. PCR for Sequencing library

- Choose the optimal PCR cycle # based on Bioanalyzer QC and amplify
- Make a PCR mix consisting of:

PCR Mix	1 reaction
water	30.8 $\mu\text{L}$
10x AccuPrime PCR buffer 1	5 $\mu\text{L}$
Primer 1 (2P_univP5, 12.5 $\mu\text{M}$ )	2 $\mu\text{L}$
Primer 2 (2P_barcode, 12.5 $\mu\text{M}$ )	2 $\mu\text{L}$
cDNA	10 $\mu\text{L}$
AccuPrime HiFi Taq (5U/ $\mu\text{L}$ )	0.2 $\mu\text{L}$
<b>Total</b>	50 $\mu\text{L}$

- Mix well and aliquot 8ul into 6 wells of a 384 well plate (amplification is more robust in smaller volumes)
- Place each place in a thermal cycler with the following conditions  
Cycling conditions:  
 1 cycle: 98°C, 3min  
 cycles: # from test PCR  
     98°C, 30sec  
     55°C, 30sec  
     65°C, 60sec  
 1cycle: 65°C, 10min  
 4°C, hold

#### 16. SPRI clean up (1.5x):

- Pool PCR reaction and increase to 2x reaction volume with sterile water
- Add 1.5x AMPure beads, and mix up/down 15x
- Incubate at RT for 15min
- Place on magnet for 5min or until solution is clear
- Aspirate out and discard clear solution
- Add 200  $\mu\text{L}$  fresh 70% EtOH without removing from magnet and incubate for 30 sec
- Aspirate out and discard the clear solution
- Add 200  $\mu\text{L}$  fresh 70% EtOH without removing from magnet and incubate for 30 sec
- Aspirate out and discard the clear solution and let air dry for 10min
- Elute off beads with 25  $\mu\text{L}$  sterile water

#### 17. Final SPRI clean up (0.7x):

- Add 0.7x AMPure beads (17.5uL), and mix up/down 15x
- Incubate at RT for 15min



- Place on magnet for 5min or until solution is clear
- Aspirate out and discard clear solution
- Add 200  $\mu$ L fresh 70% EtOH without removing from magnet and incubate for 30 sec
- Aspirate out and discard the clear solution
- Add 200  $\mu$ L fresh 70% EtOH without removing from magnet and incubate for 30 sec
- Aspirate out and discard the clear solution and let air dry for 10min
- Elute off beads with 25  $\mu$ L 1x low TE (10 mM Tris, 0.1M EDTA)

**Ordering Barcoded Oligos: (regular desalting is sufficient)**

- The barcoded oligos and 3Tr3 second adapter oligo require a 5' phosphate and 3' blocking group (either a 3'-spacer C3 or ddC)  
e.g., 3Tr3: 5'-/5Phos/AGA TCG GAA GAG CAC ACG TCT G/3SpC3/-3'
- No special modifications for the AR2 oligo (regular desalting):
  - 5' TAC ACG ACG CTC TTC CGA T 3'

**Set of 32 DNA barcoded oligos:**

barcode sequence (5'-3')	oligo	barcode read
ACATTATT	A <b>ACATTATT</b> AGATCGGAAGAGCGTCGTGTA	AATAATGT
AAGTGTTG	A <b>AAGTGTTG</b> AGATCGGAAGAGCGTCGTGTA	CAACACTT
AGAATTAT	A <b>AGAATTAT</b> AGATCGGAAGAGCGTCGTGTA	ATAATTCT
ATATGGAC	A <b>ATATGGAC</b> AGATCGGAAGAGCGTCGTGTA	GTCCATAT
ATCACTTG	A <b>ATCACTTG</b> AGATCGGAAGAGCGTCGTGTA	CAAGTGAT
CCAAGTCG	A <b>CCAAGTCG</b> AGATCGGAAGAGCGTCGTGTA	CGACTTGG
CAACTCGC	A <b>CAACTCGC</b> AGATCGGAAGAGCGTCGTGTA	GCGAGTTG
CCCGTCTT	A <b>CCCGTCTT</b> AGATCGGAAGAGCGTCGTGTA	AAGACGGG
CCCTACAG	A <b>CCCTACAG</b> AGATCGGAAGAGCGTCGTGTA	CTGTAGGG
CCCTCGGC	A <b>CCCTCGGC</b> AGATCGGAAGAGCGTCGTGTA	GCCGAGGG
CCGGTACC	A <b>CCGGTACC</b> AGATCGGAAGAGCGTCGTGTA	GGTACCGG
CGGAGGGC	A <b>CGGAGGGC</b> AGATCGGAAGAGCGTCGTGTA	GCCCTCCG
CTCGGTAC	A <b>CTCGGTAC</b> AGATCGGAAGAGCGTCGTGTA	GTACCGAG
CGGCACTT	A <b>CGGCACTT</b> AGATCGGAAGAGCGTCGTGTA	AAGTGCCG
CTCTAACT	A <b>CTCTAACT</b> AGATCGGAAGAGCGTCGTGTA	AGTTAGAG
CTGGATCG	A <b>CTGGATCG</b> AGATCGGAAGAGCGTCGTGTA	CGATCCAG
GCAGCCAC	A <b>GCAGCCAC</b> AGATCGGAAGAGCGTCGTGTA	GTGGCTGC
GAGATTGT	A <b>GAGATTGT</b> AGATCGGAAGAGCGTCGTGTA	ACAATCTC
GAGCCATC	A <b>GAGCCATC</b> AGATCGGAAGAGCGTCGTGTA	GATGGCTC
GTA ACTGC	A <b>GTA ACTGC</b> AGATCGGAAGAGCGTCGTGTA	GCAGTTAC
GGCCCAAG	A <b>GGCCCAAG</b> AGATCGGAAGAGCGTCGTGTA	CTTGGGCC
GTCTGGCG	A <b>GTCTGGCG</b> AGATCGGAAGAGCGTCGTGTA	CGCCAGAC
GGTCCTCT	A <b>GGTCCTCT</b> AGATCGGAAGAGCGTCGTGTA	AGAGGACC

TCATCGTG	<b>ATCATCGTG</b> AGATCGGAAGAGCGTCGTGTA	<b>CACGATGA</b>
TACAACAT	<b>ATACAACAT</b> AGATCGGAAGAGCGTCGTGTA	<b>ATGTTGTA</b>
TCCCGCGG	<b>ATCCCGCGG</b> AGATCGGAAGAGCGTCGTGTA	<b>CCGCGGGA</b>
TACAGATG	<b>ATACAGATG</b> AGATCGGAAGAGCGTCGTGTA	<b>CATCTGTA</b>
TACCGGCC	<b>ATACCGGCC</b> AGATCGGAAGAGCGTCGTGTA	<b>GGCCGGTA</b>
TTACCACG	<b>ATTACCACG</b> AGATCGGAAGAGCGTCGTGTA	<b>CGTGGTAA</b>
TGAACCAG	<b>ATGAACCAG</b> AGATCGGAAGAGCGTCGTGTA	<b>CTGGTTCA</b>
TGGGAGAC	<b>ATGGGAGAC</b> AGATCGGAAGAGCGTCGTGTA	<b>GTCTCCCA</b>
TTTCTAAC	<b>ATTTCTAAC</b> AGATCGGAAGAGCGTCGTGTA	<b>GTTAGAAA</b>

Set of 54 RNA barcoded oligos:

Barcode	Sequence	6 base Barcode Read + T
rUrU rGrCrU rU	/5Phos/rArUrU <b>rGrCrU</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	AAGCAAT
rUrG rArArU rU	/5Phos/rArUrG <b>rArArU</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	AATTCAT
rArC rUrUrG rU	/5Phos/rArArC <b>rUrUrG</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	ACAAGTT
rGrG rCrUrG rU	/5Phos/rArGrG <b>rCrUrG</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	ACAGCCT
rUrC rArGrG rU	/5Phos/rArUrC <b>rArGrG</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	ACCTGAT
rArU rUrArG rU	/5Phos/rArArU <b>rUrArG</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	ACTAATT
rUrU rGrArG rU	/5Phos/rArUrU <b>rGrArG</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	ACTCAAT
rUrG rGrUrC rU	/5Phos/rArUrG <b>rGrUrC</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	AGACCAT
rGrU rUrUrA rU	/5Phos/rArGrU <b>rUrUrA</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	ATAAACT
rArU rGrUrA rU	/5Phos/rArArU <b>rGrUrA</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	ATACATT
rArU rUrGrA rU	/5Phos/rArArU <b>rUrGrA</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	ATCAATT
rGrU rGrGrA rU	/5Phos/rArGrU <b>rGrGrA</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	ATCCACT
rGrA rArGrA rU	/5Phos/rArGrA <b>rArGrA</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	ATCTTCT
rGrG rUrCrA rU	/5Phos/rArGrG <b>rUrCrA</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	ATGACCT
rGrA rGrCrA rU	/5Phos/rArGrA <b>rGrCrA</b> rUrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	ATGCTCT
rCrA rGrUrU rG	/5Phos/rArCrA <b>rGrUrU</b> rGrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	CAACTGT
rGrG rUrGrU rG	/5Phos/rArGrG <b>rUrGrU</b> rGrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	CACACCT
rGrA rGrGrU rG	/5Phos/rArGrA <b>rGrGrU</b> rGrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	CACCTCT

rUrU rCrGrU rG	/5Phos/rArUrU rCrGrU rGrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	CACGAAT
rUrC rArGrU rG	/5Phos/rArUrC rArGrU rGrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	CACTGAT
rArA rGrCrU rG	/5Phos/rArArA rGrCrU rGrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	CAGCTTT
rGrA rCrUrG rG	/5Phos/rArGrA rCrUrG rGrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	CCAGTCT
rArU rCrGrG rG	/5Phos/rArArU rCrGrG rGrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	CCCGATT
rArG rUrCrG rG	/5Phos/rArArG rUrCrG rGrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	CCGACTT
rUrA rUrArG rG	/5Phos/rArUrA rUrArG rGrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	CCTATAT
rCrG rGrArG rG	/5Phos/rArCrG rGrArG rGrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	CCTCCGT
rCrU rArArG rG	/5Phos/rArCrU rArArG rGrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	CCTTAGT
rGrU rGrGrU rC	/5Phos/rArGrU rGrGrU rCrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	GACCACT
rGrG rGrArU rC	/5Phos/rArGrG rGrArU rCrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	GATCCCT
rGrA rGrUrG rC	/5Phos/rArGrA rGrUrG rCrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	GCACTCT
rArU rGrGrG rC	/5Phos/rArArU rGrGrG rCrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	GCCCATT
rGrU rArUrC rC	/5Phos/rArGrU rArUrC rCrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	GGATACT
rGrU rUrGrC rC	/5Phos/rArGrU rUrGrC rCrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	GGCAACT
rGrA rGrGrC rC	/5Phos/rArGrA rGrGrC rCrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	GGCCTCT
rCrA rArGrC rC	/5Phos/rArCrA rArGrC rCrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	GGCTTGT
rGrG rGrUrA rC	/5Phos/rArGrG rGrUrA rCrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	GTACCCT
rGrC rArArA rC	/5Phos/rArGrC rArArA rCrArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	GTTTGCT
rUrG rArUrU rA	/5Phos/rArUrG rArUrU rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TAATCAT
rArC rUrGrU rA	/5Phos/rArArC rUrGrU rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TACAGTT
rUrA rUrGrU rA	/5Phos/rArUrA rUrGrU rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TACATAT
rCrA rGrGrU rA	/5Phos/rArCrA rGrGrU rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TACCTGT
rGrU rCrGrU rA	/5Phos/rArGrU rCrGrU rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TACGACT
rArG rUrArU rA	/5Phos/rArArG rUrArU rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TATACTT
rGrU rGrArU rA	/5Phos/rArGrU rGrArU rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TATCACT
rUrU rUrUrG rA	/5Phos/rArUrU rUrUrG rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TCAAAAT

rGrA rUrGrG rA	/5Phos/rArGrA rUrGrG rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TCCATCT
rArC rArGrG rA	/5Phos/rArArC rArGrG rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TCCTGTT
rGrG rUrUrC rA	/5Phos/rArGrG rUrUrC rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TGAACCT
rArU rGrUrC rA	/5Phos/rArArU rGrUrC rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TGACATT
rArU rArGrC rA	/5Phos/rArArU rArGrC rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TGCTATT
rGrG rArCrC rA	/5Phos/rArGrG rArCrC rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TGGTCCT
rGrU rArArC rA	/5Phos/rArGrU rArArC rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TGTTACT
rCrG rGrGrA rA	/5Phos/rArCrG rGrGrA rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TTCCCGT
rGrC rGrGrA rA	/5Phos/rArGrC rGrGrA rArArG rArUrC rGrGrA rArGrA rGrCrG rUrCrG rUrGrU rArG/3SpC3/	TTCCGCT

Library Enrichment Primers (regular desalting):

Enr\_P5\_RNAtag: 5'

AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCT 3'

P7 Enrichment Primer Sequence (5' --> 3') with barcode	barcodes (BC)	BC READ (reverse complement)
CAAGCAGAAGACGGCATAACGAGATTCGTGTGCGTGACTGGAGTTCAGACGTTGCTCTTCCGATCT	TCGTGTGC	GCACACGA
CAAGCAGAAGACGGCATAACGAGATTCGCCAGAGTGACTGGAGTTCAGACGTTGCTCTTCCGATCT	TCGCCAGA	TCTGGCGA
CAAGCAGAAGACGGCATAACGAGATTCGCATATGGTGACTGGAGTTCAGACGTTGCTCTTCCGATCT	TCGCATATG	CATAGCGA
CAAGCAGAAGACGGCATAACGAGATGGCTCCTGGTGACTGGAGTTCAGACGTTGCTCTTCCGATCT	GGCTCCTG	CAGGAGCC
CAAGCAGAAGACGGCATAACGAGATATCCGACAGTGACTGGAGTTCAGACGTTGCTCTTCCGATCT	ATCCGACA	TGTCGGAT
CAAGCAGAAGACGGCATAACGAGATAACATAAATGTGACTGGAGTTCAGACGTTGCTCTTCCGATCT	AACATAAT	ATTATGTT
CAAGCAGAAGACGGCATAACGAGATATGGTAGGGTGACTGGAGTTCAGACGTTGCTCTTCCGATCT	ATGGTAGG	CCTACCAT
CAAGCAGAAGACGGCATAACGAGATGCTAAGTAGTGACTGGAGTTCAGACGTTGCTCTTCCGATCT	GCTAAGTA	TACTTAGC

## APPENDIX:

### A. Modifications of RNAtag-Seq protocol used in generation of mouse data

10ug of total RNA from mouse ES cells was fragmented using Ambion Fragmentation reagent (Invitrogen) for 2 min and 30 sec at 70°C and placed on ice. Fragmented RNA was repaired for 35min at 37°C in 400ul of FNK buffer (50mM Tris, 50mM KCl, 5mM MgCl<sub>2</sub>, 5.5mM DTT and 0.01% Triton X-100) containing 8ul of Murine RNase inhibitor (NEB), 18ul of T4 PNK (NEB) and 18ul of FastAP (Fermentas). RNA was cleaned using the Zymo RNA concentrator-5 kit (Zymo Research). 20ng of repaired RNA was ligated to 20 pmoles unique RNA barcodes (5'-P and 3'-SpC3 modifications) in separate wells. After ligation all RNA samples were pooled and cleaned using the Zymo RNA concentrator-5 kit (Zymo Research). cDNA was synthesized using

10 pmoles of AR23s primer from IDT: 5' – /5SpC3/CCTACACGACGCTCTTCC – 3' and AffinityScript Reverse Transcriptase (Agilent) at 54°C for 50min. RT primer AR23s was digested after RT using ExoSAP-IT (Affymetrix), RNA was degraded and cDNA was cleaned using Silane beads (Invitrogen). The 3'-end of cDNA was ligated with 3Tr3 adapter overnight. cDNA was amplified using the NEBNext PCR module (NEB) and indexed primers (IDT).

**B. Quick Protocol for Ribo-Zero™ Magnetic Kit (Bacteria) – For pools containing 1-5µg of total RNA. For lower input, consult Epicentre for low input protocol.**

For experience users only! – see manual for more detailed protocol

Step	Procedure
Prepare Magnetic Beads	Add 225 µl Magnetic Beads to RNase-free tube Magnetize for 1 min at RT Wash with 225 µl RNase-Free Water Magnetize, repeat wash step Resuspend in 65 µl Resuspension Solution Optional: Add 1 µl RiboGuard RNase Inhibitor
Treat sample with rRNA Removal Solution	Mix in 40 µl total volume: See manual 1-5 µg total RNA (1-2.5 µg use <b>28 µL</b> ; >2.5-5 µg use <b>26 µl</b> ) rRNA Removal Solution (1-2.5 µg use <b>8 µL</b> ; >2.5-5 µg use <b>10 µl</b> ) 4 µl Reaction Buffer Incubate 10 min @ 68°C, then 5 min @ RT
Remove rRNA	Mix previously prepared Magnetic Beads Add RNA mixture, mix well by pipetting, vortex briefly Incubate 5 min @ RT, vortex Incubate 5 min @ 50°C Magnetize, transfer supernatant (rRNA-depleted sample) to RNase-free tube (~90 µL)
Purify rRNA-depleted sample	<b>SPRI cleanup (1.8x):</b> <ul style="list-style-type: none"> <li>• Add 160 µL of AMPure XP beads containing 85-90 µL of rRNA-depleted RNA. Pipette entire volume up/down 10x</li> <li>• Incubate at RT for 15min</li> <li>• Place tubes on magnet for at least 5min (until liquid appears clear)</li> <li>• Remove and discard supernatant. Take care not to disturb the beads. Keep tubes on magnet</li> <li>• Add 200 µL freshly prepared 80% EtOH without disturbing beads, and incubate at RT for 30sec</li> <li>• Remove and discard all of the supernatant from tube, take care not to disturb the beads</li> <li>• Repeat EtOH wash (total of 2, 80% EtOH washes)</li> <li>• Air dry on magnet at RT for 15min</li> <li>• Elute off beads by mixing <b>14µL</b> of RNase/DNase free water</li> </ul>